

GM12

Gingiva Mask

Gingiva Mask GM12 is an excellent choice for replicating gingiva contours for digital implant cases. GM12 is a flexible material that can assist in planning the emergence profile for digital implant restorations.

→ Flexible, gingiva color



■ Physical Properties - Liquid Status

Color	Transparent Pink
Viscosity @25°C cps	700
Wavelength	405nm

■ Mechanical Properties - Post Cured Parts

Property	Value	Test Method
Hardness Shore A	55	ISO 868
Elongation at break %	110	ISO 527

Post Processing	Value
Cleaning time with 95% ethanol	Basic Ultrasonic: 60s + Fine Ultrasonic: 10s
Recommended curing time	FabCure 1/1S: 3 mins, FabCure: 2 mins